

## *Op-perfect excellent LATEX*

### *Op-perfect excellent*

*Latex, powder-free, examination glove optimal for tactility, ease of donning and protection*

*This non-sterile glove has a chlorinated inner surface designed to ensure easy and fast donning. Its fully textured surface provides excellent tactile sensitivity, making it ideal for patient handling and double gloving. OP-PERFECT natural rubber latex formulation provides a comfortable, well-fitting glove while extensive leaching phases and post-production washing cycles ensure a very low residual protein content. CE Marked as a Class I Medical Device, ensuring more consistent quality and providing protection for a wide range of applications.*

#### *KEY FEATURES AND BENEFITS:*

- Unique latex formulation offers comfortable, precise handling with excellent sensitivity*
- Chlorinated inner coating to ensure easy donning*
- Secure barrier protection with a high degree of tear resistance*

#### *Caution:*

- This product contains natural rubber latex which may cause allergic reactions including anaphylactic responses.*

#### *Recommended For*

- Enhanced tactile sensitivity and dexterity, especially for delicate procedures*
- Brief, low-risk examination procedures*

#### *OP-PERFECT*

*Every day, Everywhere*

*Protection*

*OP-PERFECT excellent medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include HARIR ergonomic technologies that ensure comfort during extended procedures, wet grip innovations*

## *Op-perfect excellent LATEX*

*that prevent harmful mishandling, and micro-thin films for maximum sensitivity  
Caregivers who are fully prepared and comfortably protected by versatile op-  
perfect styles have the peace of mind to focus on what matters most: the health and  
protection of themselves and their patients.*

<b>PRODUCT INFORMATION</b>			
<i>Material</i>	<i>Natural Rubber Latex</i>	<i>Not Made from Natural Rubber Latex</i>	<i>No</i>
<i>Color</i>	<i>White</i>	<i>Cuff Length</i>	<i>Standard</i>
<i>Powder Content</i>	<i>Powder-Free</i>	<i>External Glove Surface</i>	<i>micro-Textured</i>
<i>Freedom from Holes (Inspection level I)</i>	<i>1.5 AQL</i>	<i>Palm Thickness (mm)</i>	<i>0.11-012</i>
<i>Finger Thickness (mm)</i>	<i>0.12</i>	<i>Available Sizes</i>	<i>XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)</i>
<i>Tested For Use with Chemotherapy Drugs</i>	<i>No</i>	<i>Sterile</i>	<i>No</i>
<i>Antistatic</i>	<i>Not Tested</i>	<i>Glove Length (mm)</i>	<i>258-270</i>
<i>Product Segmentation</i>	<i>Essential Series</i>		

<b>ORDERING INFORMATION</b>					
<i>Size</i>	<i>XS (5.5 - 6)</i>	<i>S (6.5 - 7)</i>	<i>M (7.5 - 8)</i>	<i>L (8.5 - 9)</i>	<i>XL (9.5 - 10)</i>
<i>Product Code</i>	<i>731011133001</i>	<i>731011133002</i>	<i>731011133003</i>	<i>731011133004</i>	<i>731011133005</i>

<b>PACKAGING AND STORAGE</b>	
<i>Packaging</i>	<i>100 gloves per dispenser 10 dispensers per case 1000 gloves per case *SIZE XL</i>

## *Op-perfect excellent LATEX*

<i>Shelf Life</i>	<i>5 Years</i>
<i>Storage Instructions</i>	<i>Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.</i>

### **STANDARDS AND CERTIFICATIONS**

*EN 374-1:2016, EN 21420:2020, , EN 455 1-4, EN ISO 374-5:2016, ISO 10993-10:2021, ISO 13485, ISO 9001, ISO 10002 :2018, ISO 10004:2018, HALAL, FDA QSR 21 CFR 820 (c GMP), GDP Good Distribution Practices*

*\*Note: If this size is requested, the number of gloves will be based on their weight and length.*