

Op-perfect premium LATEX

Op-perfect premium

Latex, low-powder, examination glove with premium grip properties

This non-sterile glove has excellent grip levels making it ideal for healthcare providers who handle instrumentation, test tubes or infusion equipment. OP-PERFECT natural rubber latex formulation provides a comfortable, well-fitting glove while the smooth powdered outer surface facilitates ease of donning.

KEY FEATURES AND BENEFITS:

- Textured & powdered outer surface provides an excellent grip in wet and dry conditions*
- Smooth powdered outer surface for easy donning*
- Low allergenicity latex reduces risk of contact dermatitis due to absence of thiurams*

Caution:

- This product contains natural rubber latex which may cause allergic reactions including anaphylactic responses.*
- surface powder should be aseptically removed prior to use*

Recommended For

- Brief, low-risk examination procedures*

OP-PERFECT

Every day, Everywhere

Protection

OP-PERFECT premium medical

OP-PERFECT premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Harir ergonomic technologies that ensure comfort during extended procedures, powdered for ease donning, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are

Op-perfect premium LATEX

fully prepared and comfortably protected by versatile MICROTOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

PRODUCT INFORMATION

<i>Material</i>	<i>Natural Rubber Latex</i>	<i>Not Made from Natural Rubber Latex</i>	<i>No</i>
<i>Color</i>	<i>White</i>	<i>Cuff Length</i>	<i>Standard</i>
<i>Powder Content</i>	<i>Low powder</i>	<i>External Glove Surface</i>	<i>micro-Textured</i>
<i>Freedom from Holes (Inspection level I)</i>	<i>1.5 AQL</i>	<i>Palm Thickness (mm)</i>	<i>0.11-012</i>
<i>Finger Thickness (mm)</i>	<i>0.11-012</i>	<i>Available Sizes</i>	<i>XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)</i>
<i>Tested For Use with Chemotherapy Drugs</i>	<i>No</i>	<i>Sterile</i>	<i>No</i>
<i>Antistatic</i>	<i>Not Tested</i>	<i>Glove Length (mm)</i>	<i>255-267</i>
<i>Product Segmentation</i>	<i>Essential Series</i>		

ORDERING INFORMATION

<i>Size</i>	<i>XS (5.5 - 6)</i>	<i>S (6.5 - 7)</i>	<i>M (7.5 - 8)</i>	<i>L (8.5 - 9)</i>	<i>XL (9.5 - 10)</i>
<i>Product Code</i>	<i>731011133001</i>	<i>731011133002</i>	<i>731011133003</i>	<i>731011133004</i>	<i>731011133005</i>

PACKAGING AND STORAGE

<i>Packaging</i>	<i>100 gloves per dispenser 10 dispensers per case 1000 gloves per case *SIZE XL</i>
<i>Shelf Life</i>	<i>5 Years</i>

Op-perfect premium LATEX

<i>Storage Instructions</i>	<i>Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.</i>
---------------------------------	---

STANDARDS AND CERTIFICATIONS

EN 374-1:2016, EN 21420:2020, EN 455 1-4, EN ISO 374-5:2016, ISO 10993-10:2021, ISO 13485, ISO 9001, ISO 10002 :2018, ISO 10004:2018, HALAL, FDA QSR 21 CFR 820 (c GMP), GDP Good Distribution Practices

**Note: If this size is requested, the number of gloves will be based on their weight and length.*