

Packing plan for NITEX

All issues - warning notices, signs are used, identification plate, storage, transportation, and supply conditions in the plan package are available.

Product Name: nitrile examination gloves

Brand: Nitex

1) Packing items:

The initial packaging of the examination gloves is in a dispenser.

1.1) Packaging information:

NAME	<i>Nitex 100 *20 PCS Size L</i>	<i>Nitex 100*20 PCS Size M</i>	<i>Nitex 100*20 PCS Size S</i>
<i>Dispenser WEIGHT (GROSS) (100 pcs)</i>	<i>500 gr</i>	<i>465 gr</i>	<i>405</i>
<i>Case WEIGHT (GROSS) (2000 pcs)</i>	<i>10800 gr</i>	<i>9600 gr</i>	<i>8800 gr</i>
<i>Case DIMENSION (LENGTH)</i>	<i>450 mm</i>	<i>450 mm</i>	<i>450 mm</i>
<i>Case DIMENSION (WIDTH)</i>	<i>260 mm</i>	<i>260 mm</i>	<i>260 mm</i>
<i>Case DIMENSION (HEIGHT)</i>	<i>260 mm</i>	<i>260 mm</i>	<i>260 mm</i>
<i>All weights are approximate. make sure to get detailed information from the seller before shipment</i>			

1.2) Number in dispenser: 100 pcs

2) product Information:

-nitrile examination glove

-Powder free

-Excellent finishing

-medical grade


-AQL 1.5



- Type of product: **Examination Gloves** • Material of product: **Nitrile**
- Specification: **Non-Sterile, Powder Free, Easy wear, Ambidextrous, Micro textured, Single use only**
- Number in a box: **100 Pieces by weight**
- Storage: **Store at room temperature, avoid excessive heat (5 - 40° C), Pack should be shield from, sunlight, fluorescent lighting, X-ray, moisture and ozone**
- Germany Identifier: **SONNHEIM Medikal (Germany Reg. No. 39634985.4)**
- Manufacturer name & address: **Iran Harir Glove Co. Urmia Factory: Phase #2 of Industrial city, Urmia - I.R.I**
- **EC REP TAI Trading Agency GmbH**
Schüttelstraße 55/11, 1020 Wien-AUSTRIA
Tel.: +43 (0) 732 77 33 20 / Fax: +43 (0) 732 77 33 20-15
E-Mail: info@tai-austria.com
- Customer Service (Medical Dep.) Tel: + 98 (21) 88 91 64 36 Fax: + 98 (21) 88 91 64 37
www.glove.ir / medical@glove.ir Zip Code: 1415846351 Made in I.R.I

2.2) CAUTION: -

2.3) symbol used:



3) information for transportation:

Product name	Capacity on pallet (120 * 100 cm)	Capacity in transit trailer, without pallet (approx.)	Capacity in transit trailer, with pallet (approx.)	Capacity in 40-foot container, without pallet (approx.)	Capacity in: 40-foot container, with pallet (approx.)
nitex Examination glove (20*100)	72	XS/S/M=2600 L=2450	1872	2000	1440

Size M

Excellent Finishing Powder Free Nitrile Examination Glove



Excellent Finishing Powder Free Nitrile Examination Glove
Size M

NITEX

Excellent Finishing Powder Free Nitrile Examination Glove

Latex Free فاقد لاتکس

100 Pieces by Weight

Size M

CE, ISO 9001, ISO 13485, ISO 14001, ISO 45001, ISO 22716, ISO 22717

حلال HALAL

Size M
Excellent Finishing Powder Free Nitrile Examination Glove

NITEX

Excellent Finishing Powder Free Nitrile Examination Glove

Size M

NITEX

Iran Harir Glove Co. Since 1977

Multi Purpose Disposable Glove

Excellent Finishing Powder Free Nitrile Examination Glove

Size M

Medical Grade AQL 1.5

- XS
- S
- M
- L
- XL
- XXL

- XS
- S
- M
- L
- XL
- XXL

• Type of product: Examination Gloves • Material of product: Nitrile
 • Specification: Non-Sterile, Powder Free, Easy wear, Antidistress, Micro textured, Single use only
 • Number in a box: 100 Pieces by weight
 • Storage: Store at room temperature, avoid excessive heat (5 - 40° C).
 Pack should be shield from, sunlight, fluorescent lighting, X-ray, moisture and ozone
 • Germany Identifier: SCANNHEIM Medikal (Germany Reg. No. 39634985.4)
 • Manufacturer name & address: Iran Harir Glove Co.
 Urmia Factory, Phase #2 of Industrial city, Urmia - I.R.I.
TAI Trading Agency GmbH
 Schallstrasse 55/1, 1020 WIEN-AUSTRIA
 Tel.: +43 (0) 732 77 33 20 / Fax: +43 (0) 732 77 33 20-15
 E-Mail: info@tai-austria.com
 • Customer Service (Medical Dep.) Tel: + 98 (21) 88 91 64 36 Fax: + 98 (21) 88 91 64 37
 www.gloves.ir / medical@gloves.ir Zip Code: 1415646381 Made in I.R.I.



CE
 ISO 9001
 ISO 13485
 ISO 10002 & 10004
 FDA CFR 21 CFR 820 (cGMP)
 GDP Good Distribution Practices
 CE (Directive 93/42/EEC, Class I)

شماره پیامک استاندارد ۱۰۰۰۱۵۱۷

توزیع کننده ارومیه: شهر صنعتی، فاز ۲ - پ. ۶۱۷۶۶۱۹۳
 توزیع کننده ارومیه: شهر صنعتی، فاز ۲، خیابان حافظ هشتاد، پ. ۱۰۰۰۱۵۱۷
 پ. ۶۶۶۶۶۶۶۶ • بانک پارس کوثر
 ساخت ایران

عزم همه مایه و باطن برای انجام علوم پزشکی است. هرگز از این دستکش برای کارهای غیر پزشکی استفاده نکنید.
 هرگز از دستکش برای کارهای غیر پزشکی استفاده نکنید.