

Packing plan for OP-PERFECT premium

All issues- warning notices, signs are used, identification plate, storage, transportation, and supply conditions in the plan package are available.

Product Name: low powder latex examination gloves

Brand: OP-PERFECT premium

1) Packing items:

The initial packaging of the examination gloves is in a Dispenser.

1.1) Packaging information:

NAME	LATEX 100 PCS Size L	LATEX 100 PCS Size M	LATEX 100 PCS Size S
Dispenser WEIGHT (GROSS) (100 pcs)	646 gr	574 gr	550 gr
Case WEIGHT (GROSS) (1000 pcs)	7000 gr	6330 gr	5980 gr
Case DIMENSION (LENGTH)	405mm	405mm	405mm
Case DIMENSION (WIDTH)	255 mm	255 mm	255 mm
Case DIMENSION (HEIGHT)	235mm	235mm	235mm
All weights are approximate. make sure to get detailed information from the seller before shipment.			

1.2) Number in Dispenser: 100 pcs

2) product Information:

- latex examination glove*
- low powder*
- Excellent finishing*
- medical grade*
- AQL 1.5*



- Type of product: **Examination Gloves** ● Material of product: **Latex (NR)**
- Specification: **Non-Sterile, Easy wear, Polymer coated low powder, Ambidextrous, Micro textured, Single use only**
- Number in a box: **100 Pieces by weight**
- △ Cautions: **This product contains natural rubber latex which may cause allergic reactions in some individuals**
- Surface powder should be removed after donning
- Storage: **Store at room temperature, avoid excessive heat (5 - 40° C), Pack should be shield from sunlight, fluorescent lighting, X-ray, moisture and ozone**
- Germany Identifier: **SONNHEIM Medikal (Germany Reg. No. 39634985.4)**
- Manufacturer name & address: **Iran Harir Glove Co. / Urmia Factory: Phase #2 of Industrial city, Urmia - I.R.I**
- EC REP **TAI Trading Agency GmbH / Sch ü ttelstra ß e 55/11, 1020 Wien-AUSTRIA / Tel.: +43 (0) 732 77 33 20**
- Fax: **+43 (0) 732 77 33 20-15 / E.Mail: info@tai-austria.com**
- Customer Service (Medical Dep.) Tel: **+ 98 (21) 88 91 64 36** Fax: **+ 98 (21) 88 91 64 37**
- www.glove.ir / medical@glove.ir** Zip Code: **1415846351** Made in I.R.I

2.2) CAUTION:

△ Cautions: **This product contains natural rubber latex which may cause allergic reactions in some individuals**

2.3) symbol used:



3) information for transport:

Product name	Capacity on pallet (110 * 100 cm)	Capacity in transit trailer, without pallet (approx.)	Capacity in transit trailer, with pallet (approx.)	Capacity in 40-foot container, without pallet (approx.)	Capacity in: 40-foot container, with pallet (approx.)
low powder Examination glove	90	3260	2340	2860	1980

OP-PERFECT Low Powder Examination Glove
PREMIUM FINISHING
Size **M**

Latex
Low Powder

CE 100 Pieces by Weight

ISO 9001
ISO 13485
ISO 15002 & 19004
FDA QSR 21 CFR 820 (cGMP)
GMP Good Distribution Practices
CE (Directive 93/42/EEC, Class I)

Iran Harir Glove Co.
Since 1977

OP-PERFECT Low Powder Examination Glove
PREMIUM FINISHING
Size **M**

Medical Grade
AQL 1.5

6 260045 228095

• Type of product: Examination Gloves • Material of product: Latex (NR)
• Specification: Non-Sterile, Easy wear, Polymer coated low powder, Ambidextrous, Micro textured, Single use only
• Number in a box: 100 Pieces by weight
⚠️ Cautions: This product contains natural rubber latex which may cause allergic reactions in some individuals
• Surface powder should be removed after donning
• Storage: Store at room temperature, avoid excessive heat (5 - 40° C), Pack should be shield from sunlight, fluorescent lighting, X-ray, moisture and ozone
• Germany Identifier: SONNHEIM Medikal (Germany Reg. No. 39634985 4)
• Manufacturer name & address: Iran Harir Glove Co. / Urmia Factory: Phase #2 of Industrial city, Urmia - I.R.I.
• Trading Agency GmbH / Sch. u. Telstra • 5511, 1100 Wien-AUSTRIA / TEL: +43 (0) 732 77 33 20
FAX: +43 (0) 732 77 33 20-15 / E-Mail: info@tai-austria.com
• Customer Service (Medical Dep.) Tel: +98 (21) 88 91 64 36 Fax: +98 (21) 88 91 64 37
www.glove.ir / medical@glove.ir Zip Code: 1415846351 Made in I.R.I.