

# *NITEX examination nitril gloves*

## *NITEX*

*Nitrile Examination Glove with Textured Fingertips  
Offers comfortable and reliable protection*

- *Designed with a lightweight nitrile formulation for enhanced comfort and tactile sensitivity*
- *Textured on the fingertips to promote sure handling of materials and instruments*
- *Packaged with 100 gloves per dispenser*
- *Offers reliable protection against hazardous chemicals, including chemotherapy drugs*

## *Industries*

- *Healthcare*
- *Emergency Medical Services*
- *Life Sciences*
- *Dental sciences/dental industries*

## *Recommended For*

- *Standard patient examinations*
- *General dentistry*
- *Emergency services*
- *Laboratory analysis*
- *Protection from Type I latex allergy in HCW's or patients*

# *NITEX examination nitril gloves*

## *TECHNICAL DATA SHEET*

<i>PRODUCT INFORMATION</i>			
<i>Material</i>	<i>Nitrile</i>	<i>Manufacturing/QMS Audit Standards</i>	<i>ISO 13485, ISO 9001</i>
<i>Color</i>	<i>Blue (Multicolor on request)</i>	<i>Cuff</i>	<i>Beaded</i>
<i>Glove Design</i>	<i>Chlorinated, Textured Fingers</i>	<i>Regulatory/Standards Compliance</i>	<i>ASTM D6319, EN 16523-1, EN 21420:2020, EN 455 1-4, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 13485, ISO 9001</i>
<i>Packaging</i>	<i>100 gloves per dispenser; 20 dispensers/2000 gloves by weight</i>	<i>Storage</i>	<i>Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition</i>
<i>User Needs Segment</i>	<i>High Touch</i>	<i>Available Sizes</i>	<i>XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)</i>
<i>Glove Length (mm)</i>	<i>255-267</i>	<i>Sterile</i>	<i>No</i>
<i>Antistatic</i>	<i>Not Tested</i>		

<i>PHYSICAL PROPERTIES</i>		
	<i>Typical Values</i>	<i>Testing Method</i>
<i>Length (mm)</i>	<i>255-267</i>	<i>EN 455-2</i>
<i>Freedom from Holes (Inspection level I)</i>	<i>1.5 AQL</i>	<i>EN 455-1</i>
<i>Palm Thickness (mm)</i>	<i>0.08-0.09</i>	<i>EN 455-2</i>
<i>Finger Thickness (mm)</i>	<i>0.11-0.12</i>	<i>EN 455-2</i>

## *NITEX examination nitril gloves*

### *ORDERING INFORMATION*

	<i>BEFORE AGING</i>	<i>AFTER AGING</i>	
<i>Elongation</i>	650	630	
<i>Force at break (N)</i>	$\geq 7$	$\geq 6$	

### *ORDERING INFORMATION*

<i>Size</i>	<i>XS (5.5 - 6)</i>	<i>S (6.5 - 7)</i>	<i>M (7.5 - 8)</i>	<i>L (8.5 - 9)</i>	<i>XL (9.5 - 10)</i>
<i>Product Code</i>	731032933001	731032933002	731032933003	731032933004	731032933005