

# ***OP-PERFECT vinyl gloves***

## ***OP-PERFECT vinyl***

*Vinyl, powder-free examination glove for a wide range of applications*

- *Quick and easy donning*
- *100% synthetic to reduce risk of latex related allergies or irritation*
- *Smooth surface provides a more natural feel*
- *Beaded cuff to prevent cuff roll down*

### ***Industries***

- *Healthcare*

### ***Recommended For***

- *Brief, low-risk examination procedures*
- *Protection from Type I latex allergy in HCW's or patients*

## ***TECHNICAL DATA SHEET***

<i><b>PRODUCT INFORMATION</b></i>	
<i>Material</i>	<i>Vinyl (Polyvinyl Chloride)</i>
<i>Color</i>	<i>NATURAL COLOR (Multicolor on request)</i>
<i>Glove Design</i>	<i>Polyurethane Coated, Smooth, powder free</i>
<i>Cuff</i>	<i>Beaded</i>
<i>Manufacturing/QMS Audit Standards</i>	<i>ISO 13485, ISO 9001</i>
<i>Regulatory/Standards Compliance</i>	<i>EN ISO 21420:2020 EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 9001, ISO 13485, ISO 10993-10:2021</i>
<i>Packaging</i>	<i>100 gloves per dispenser 10 dispensers per case 1000 gloves per case</i>
<i>Available sizes</i>	<i>XS (5.5 - 6), S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)</i>
<i>Anti-static</i>	<i>Not Tested</i>

## *OP-PERFECT vinyl gloves*

### *PHYSICAL PROPERTIES*

	<i>Typical Values</i>	<i>Testing Method</i>
<i>Length (mm)</i>	<i>240</i>	<i>ASTM D3767, EN 21420</i>
<i>Freedom from Holes (Inspection level I)</i>	<i>1.5 AQL</i>	<i>ASTM D5151, EN 455-1</i>
<i>Palm Thickness (mm)</i>	<i>0.08-0.09</i>	<i>ASTM D3767, EN 21420</i>
<i>Finger Thickness (mm)</i>	<i>0.11-0.12</i>	<i>ASTM D3767, EN 21420</i>

### *ORDERING INFORMATION*

	<i>BEFORE AGING</i>	<i>AFTER AGING</i>
<i>Elongation</i>	<i>350</i>	<i>350</i>
<i>Force at break (N)</i>	<i>≥ 7</i>	<i>≥ 7</i>

### *PACKAGING AND STORAGE*

<i>Packaging</i>	<i>100 gloves per dispenser 10 dispensers per case 1000 gloves per case *SIZE XL</i>
<i>Shelf Life</i>	<i>5 Years</i>
<i>Storage Instructions</i>	<i>Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.</i>

*\*Note: If this size is requested, the number of gloves will be based on their weight and length.*

### *FEATURED TECHNOLOGIES*









*NO NATURAL RUBBER*

# *OP-PERFECT vinyl gloves*

## *STANDARDS AND CERTIFICATIONS*

- 
- *EN ISO 21420:2020*

- 
-  EN 374
  -  EN 374
  -  EN ISO 374-5/Type C
  -  ISO
  -  ISO 9001
  -  ISO 13485
  - *ISO 10993-10:2002*