

## Specification for Examination Nitex gloves

<i>Color</i>	<i>Blue</i>												
<i>Shape</i>	<i>Straight Fingers</i>												
<i>Cuff</i>	<i>Beaded</i>												
<i>Powder free Nitex Examination Gloves</i>													
<i>PHYSICAL DIMENSIONS</i>													
<i>Size</i>	<i>Length (mm)</i>		<i>Palm Width (mm)</i>				<i>Thickness(mm)</i>						
	<i>Standard</i>	<i>Test</i>	<i>standard</i>	<i>Test</i>	<i>Test</i>								
					<i>Cuff</i>			<i>Palm</i>			<i>Fin</i>		
					<i>Smooth area 0.08 Min</i>			<i>Textured area 0.08 Min</i>			<i>Textured area 0.11 Min</i>		
<i>Small</i>	$\geq 230$	255-257	$80 \pm 10$	87-89	0.08-0.09			0.08-0.09			0.11-0.12		
<i>Medium</i>		261-263	$95 \pm 10$	97-99	0.08-0.09			0.08-0.09			0.11-0.12		
<i>Large</i>		265-267	$110 \pm 10$	108-110	0.08-0.09			0.08-0.09			0.11-0.12		
<i>Tensile Properties</i>		<i>Standard</i>		<i>Test</i>									
		<i>EN</i>	<i>ISO</i>	<i>1</i>	<i>2</i>	<i>3</i>	<i>4</i>	<i>5</i>	<i>6</i>	<i>7</i>	<i>8</i>	<i>9</i>	<i>10</i>
<i>Before aging</i>	<i>Max force (N)</i>	3.6	7	7.2	7.3	7.3	7.4	7.2	7.1	7.2	7.2	7.4	7.3
	<i>Elongation (%)</i>	500	650	655	658	656	653	659	655	655	658	655	658
<i>After aging</i>	<i>Max force (N)</i>	3.6	6	6.5	6.5	6.4	6.5	6.4	6.3	6.3	6.4	6.4	6.4
	<i>Elongation (%)</i>	400	500	620	620	630	620	630	620	630	620	630	620