

Specification for low powder Examination gloves Op-perfect premium

<i>Color</i>		<i>white</i>											
<i>Shape</i>		<i>Straight Fingers</i>											
<i>Cuff</i>		<i>Beaded</i>											
<i>Low powder latex Examination Gloves</i>													
<i>PHYSICAL DIMENSIONS</i>													
<i>Size</i>		<i>Length (mm)</i>		<i>Palm Width (mm)</i>				<i>Thickness(mm)</i>					
		<i>Standard</i>	<i>Test</i>	<i>standard</i>	<i>Test</i>	<i>Test</i>							
						<i>Cuff</i>		<i>Palm</i>		<i>Fin</i>			
						<i>Smooth 0.08 Min</i>		<i>Textured area 0.11 Min</i>		<i>Textured area 0.11 Min</i>			
<i>Small</i>	≥230	256-268	80±10	84-87		0.08-0.09		0.11-0.12		0.11-0.12			
<i>Medium</i>		256-268	95±10	96-99		0.08-0.09		0.11-0.12		0.11-0.12			
<i>Large</i>		256-268	110±10	106-108		0.08-0.09		0.11-0.12		0.11-0.12			
<i>Tensile Properties</i>		<i>Standard</i>		<i>Test</i>									
		<i>EN</i>	<i>ISO</i>	<i>1</i>	<i>2</i>	<i>3</i>	<i>4</i>	<i>5</i>	<i>6</i>	<i>7</i>	<i>8</i>	<i>9</i>	<i>10</i>
<i>Before aging</i>	<i>Max force (N)</i>	3.6	7	7.4	7.3	7.4	7.2	7.3	7.4	7.3	7.4	7.2	7.3
	<i>Elongation (%)</i>	500	650	680	680	690	670	670	660	670	680	670	680
<i>After aging</i>	<i>Max force (N)</i>	3.6	6	6.4	6.5	6.5	6.4	6.5	6.4	6.5	6.3	6.5	6.5
	<i>Elongation (%)</i>	400	500	630	620	620	630	640	630	640	620	640	630
<i>Protein (µgr/gr)</i>		50	50	41	43	42	42	41	40	42	43	42	43
<i>Powder (mg)</i>		180	180	60	63	63	65	63	65	65	63	65	62