

Packing plan for OP-PERFECT EXCELLENT

All issues- warning notices, signs are used, identification plate, storage, transportation, and supply conditions in the plan package are available.

Product Name: latex surgical gloves

Brand: OP-PERFECT EXCELLENT

1) Packing items:

The initial packaging of the surgical gloves is in an envelope.

1.1) Packaging information:

NAME	<i>LATEX 50 pairs Size 6.5</i>	<i>LATEX 50 pairs Size 7</i>	<i>LATEX 50 pairs Size 7.5</i>	<i>LATEX 50 pairs Size 8</i>	<i>LATEX 50 pairs Size 8.5</i>
<i>Dispenser WEIGHT (GROSS) 50 pairs</i>	<i>1740 gr</i>	<i>1780gr</i>	<i>1850 gr</i>	<i>1960 gr</i>	
<i>Case WEIGHT (GROSS) 400 pairs</i>	<i>14670 gr</i>	<i>15000gr</i>	<i>15600 gr</i>	<i>16480 gr</i>	
<i>Case DIMENSION (LENGTH)</i>	<i>620mm</i>	<i>620mm</i>	<i>620mm</i>	<i>620mm</i>	
<i>Case DIMENSION (WIDTH)</i>	<i>296 mm</i>	<i>296 mm</i>	<i>296 mm</i>	<i>296 mm</i>	
<i>Case DIMENSION (HEIGHT)</i>	<i>400mm</i>	<i>400mm</i>	<i>400mm</i>	<i>400mm</i>	

All weights are approximate. make sure to get detailed information from the seller before shipment.

1.2) Number in envelope: 1 pair

2) product Information:

-latex surgical glove

- powder free

-Excellent finishing

-medical grade

- **Type of product:** Surgical Gloves
- **Material of product:** Latex (NR)
- **Specification:** Sterile, Latex & Powder Free, Easy wear, **Polymer coated plus** Single use only • Palm textured
- The **protein** level of this product is lower than ASTM D 5712-15 Standard
- Every glove 100% is **air inflation** tested and **visually** inspected for quality and safety
- Improved **efficiency** as less glove are wasted
- Recommended for **all surgical procedures**
- **Extra Length** for more protection
- Disposable, Soft, Sensible
- **Anatomic** shape
- **Beaded** cuff
- Improved **grip** to prevent slipping **surgical instrument**
- **Sterilized by:** Ethylene oxide (**E.O**) gas
- Number in a pack: **1 Pair**
- ⚠ **Cautions:** This product contains natural rubber latex which may cause allergic reactions in some individuals
- ⊗ Sterile unless package is opened or damaged
- **Storage:** Store at room temperature, avoid excessive heat (5 - 40°C), Pack should be shield from sunlight, fluorescent lighting, x-ray, moisture and ozone
- **Identifiers:** ISO 13485 / ISO 9001 / ISO 10002 & 10004 / FDA QSR 21 CFR 820 (cGMP) GDP Good Distribution Practices / ISIRI (Cer. No. 6117041877) Health sign of Ministry of Health (Cer. No. 68937701) / SONNHEIM Medikal (Germany Reg. No. 39634985.4)
- 🏭 **Manufacturer name & address:** Iran Harir Glove Co. / Urmia Factory: Phase #2 of Industrial city, Urmia - I.R.I
- 🇪🇺 **REGP** TAI Trading Agency GmbH / Sch ü ttelstra e 55/11, 1020 Wien-AUSTRIA / Tel.: +43 (0) 732 77 33 20 Fax: +43 (0) 732 77 33 20-15 / E.Mail: info@tai-austria.com
- **Customer Service (Medical Dep.)** Tel: + 98 (21) 88 91 64 36 Fax: + 98 (21) 88 91 64 37 www.glove.ir / medical@glove.ir / **Zip Code:** 1415846351 / **Made in I.R.I**

2.2) CAUTION:

- ⚠ **Cautions:** This product contains natural rubber latex which may cause allergic reactions in some individuals
- ⊗ Sterile unless package is opened or damaged

2.3) symbol used:



1Pair
STERILE EO

Size

LOT



3) information for transportation:

<i>Product name</i>	<i>Capacity on pallet (110 * 130 cm)</i>	<i>Capacity in transit trailer, without pallet (approx.)</i>	<i>Capacity in transit trailer, with pallet (approx.)</i>	<i>Capacity in 40-foot container, without pallet (approx.)</i>	<i>Capacity in: 40-foot container, with pallet (approx.)</i>
<i>Latex surgical gloves</i>	35	1050	770	756	630



Latex, Powder Free Surgical Gloves
OP-PERFECT®
 EXCELLENT FINISHING



Iran Harir Glove Co.
 Since 1977



- **Type of product:** Surgical Gloves
- **Material of product:** Latex (NR)
- **Specification:** Sterile, Latex & Powder Free, Easy wear, Polymer coated plus Single use only • Palm textured
- The **protein** level of this product is lower than ASTM D 5712-15 Standard
- Every glove 100% is **air inflation** tested and **visually** inspected for quality and safety
- Improved **efficiency** as less glove are wasted
- Recommended for **all surgical procedures**
- **Extra Length** for more protection
- Disposable, Soft, Sensible
- **Anatomic** shape
- **Beaded** cuff
- Improved **grip** to prevent slipping **surgical instrument**
- **Sterilized by:** Ethylene oxide (E.O) gas
- Number in a pack: **1 Pair**

⚠ **Cautions:** This product contains natural rubber latex which may cause allergic reactions in some individuals
 ☉ Sterile unless package is opened or damaged

• **Storage:** Store at room temperature, avoid excessive heat (5 - 40°C). Pack should be shield from sunlight, fluorescent lighting, x-ray, moisture and ozone

• **Identifiers:** ISO 13485 / ISO 9001 / ISO 10002 & 10004 / FDA QSR 21 CFR 820 (cGMP)
 GDP Good Distribution Practices / ISIRI (Cer. No. 6117041877)

Health sign of Ministry of Health (Cer. No. 68937701) / SONNHEIM Medikal (Germany Reg. No. 39634985.4)
 • **Manufacturer name & address:** Iran Harir Glove Co. / Urmia Factory: Phase #2 of Industrial city, Urmia - I.R.I

• **Customer Service (Medical Dep.) Tel:** + 98 (21) 88 91 64 36 **Fax:** + 98 (21) 88 91 64 37
 www.glove.ir / medical@glove.ir / Zip Code: 1415846351 / **Made in I.R.I**

شماره پیامک استاندارد پروانه استاندارد ۱۰۰۰۱۵۱۷
 ساخت ایران
 • یافت ریز شد آموزش در کف دست
 حکمت هشتم، بی دست جراحی
 کارخانه قزوین: شهر صنعتی البرز، خیابان
 کارخانه ارومیه: شهر صنعتی، فاز ۲



Size

LOT



معاف از مالیات بر ارزش افزوده