

# *Op-perfect excellent SURGICAL GLOVES*

## *Op-perfect excellent*

*Latex, powder-free surgical glove with a micro-textured finish for a safe and secure grip*

*Trusted around the world to deliver the every day protection for surgeons, O.R. nurses and surgical technicians can depend on to perform their daily tasks safely and comfortably.*

### *KEY FEATURES AND BENEFITS:*

- Proprietary glove formulation delivers a soft and sensitive, yet robust surgical glove*
- Good handling and grip achieved through a unique external surface treatment*
- Micro-textured finish on the palm and fingers provides secure handling of instruments in dry and damp conditions*

*CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.*

### *Recommended For*

*Sterile procedures in the O.R.*

*Extra comfort during long, sustained procedures*

*Brief, low-risk examination procedures*

## *Op-perfect excellent SURGICAL GLOVES*

### **PRODUCT INFORMATION**

<i>Material</i>	<i>Natural Rubber Latex</i>	<i>Not Made from Natural Rubber Latex</i>	<i>No</i>
<i>Color</i>	<i>Natural color</i>	<i>Cuff Style</i>	<i>Beaded</i>
<i>Allergy Prevention</i>	<i>None</i>	<i>Sterile</i>	<i>yes</i>
<i>Powder Content</i>	<i>Powder-Free</i>	<i>outer Surface</i>	<i>micro-Textured</i>
<i>Freedom from Holes (Inspection level I)</i>	<i>0.65 AQL</i>	<i>Palm Thickness (mm)</i>	<i>0.17-018</i>
<i>Finger Thickness (mm)</i>	<i>0.23-0.24</i>	<i>Available Sizes</i>	<i>6.5-7-7.5-8-8.5</i>
<i>Inner Surface</i>	<i>Polymer Coated</i>	<i>Glove Length (mm)</i>	<i>285-305</i>
<i>Product Standards</i>	<i>EN 455 1-4, ISO 10282</i>	<i>Quality / Environmental Standards</i>	<i>EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001</i>

### **ORDERING INFORMATION**

<i>Size</i>	<i>6.5</i>	<i>7</i>	<i>7.5</i>	<i>8</i>	<i>8.5</i>
<i>Product Code</i>	<i>741031132002</i>	<i>741031132003</i>	<i>741031132004</i>	<i>741031132005</i>	<i>741031132006</i>

### **PACKAGING AND STORAGE**

<i>Packaging</i>	<i>400 pairs per case; 8 dispensers per case; 50 pairs/envelopes per dispenser</i>
<i>Shelf Life</i>	<i>5 Years</i>
<i>Storage Instructions</i>	<i>Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition</i>