

Specification for Low powder LATEX surgical gloves

<i>Color</i>	<i>Natural color</i>												
<i>Shape</i>	<i>Curve Fingers</i>												
<i>Cuff</i>	<i>Beaded</i>												
<i>Low powder latex surgical Gloves</i>													
PHYSICAL DIMENSIONS													
<i>Size</i>	<i>Length (mm)</i>			<i>Palm Width (mm)</i>				<i>Thickness(mm)</i>					
	<i>Standard</i>	<i>Test</i>	<i>Standard</i>	<i>Test</i>	<i>Standard</i>								
					<i>Cuff</i>			<i>Palm</i>			<i>Fin</i>		
					<i>Smooth 0.10 Min</i>			<i>Thextured 0.13 Min</i>			<i>Thextured 0.13 Min</i>		
					<i>Test</i>								
6.5	≥260	285	83±5	85	0.15	0.17	0.22						
7	≥270	289	89±5	92	0.14	0.16	0.22						
7.5	≥270	289	95±5	99	0.15	0.17	0.23						
8	≥270	287	102±6	106	0.14	0.16	0.22						
<i>Tensile Properties</i>		<i>Standard</i>			<i>Test</i>								
		<i>EN</i>	<i>ISO</i>	<i>1</i>	<i>2</i>	<i>3</i>	<i>4</i>	<i>5</i>	<i>6</i>	<i>7</i>	<i>8</i>	<i>9</i>	<i>10</i>
<i>Before aging</i>	<i>Max force (N)</i>	9	12.5	18.8	18.7	18.9	19	19.2	18.8	19	18.8	18.7	19
	<i>Elongation (%)</i>	-	700	1080	1080	1090	1090	1080	1090	1090	1090	1090	1080
	<i>Initial Modulus (MPa)</i>	-	2	0.48	0.50	0.48	0.49	0.48	0.48	0.48	0.48	0.49	0.48
<i>After aging</i>	<i>Max force (N)</i>	9	9.5	16.5	16.4	16.4	16.3	16.5	16.4	16.3	16.4	16.4	16.4
	<i>Elongation (%)</i>	-	550	950	960	960	960	950	950	960	960	960	950
<i>Protein (µgr/gr)</i>		50	50	34.5	35.2	34.1	35.2	33.9	34.1	35.2	35.3	34.3	36
<i>Powder (mg)</i>		>2	117	59	60	62	61	65	59	61	62	63	64