

### **Packing plan for OP-PERFECT**

All issues- warning notices, signs are used, identification plate, storage, transportation, and supply conditions in the plan package are available.

**Product Name:** vinyl examination gloves

**Brand:** OP-PERFECT

#### **1) Packing items:**

The initial packaging of the examination gloves is in a dispenser.

#### **1.1) Packaging information:**

<b>NAME</b>	<b>VINYL Size L</b>	<b>VINYL Size M</b>	<b>VINYL Size S</b>
<i>Dispenser WEIGHT (GROSS) (100 pcs)</i>	<i>594gr</i>	<i>577gr</i>	<i>554gr</i>
<i>Case WEIGHT (GROSS)(1000 pcs)</i>	<i>6400gr</i>	<i>6200gr</i>	<i>5750gr</i>
<i>Case DIMENSION (LENGTH)</i>	<i>401mm</i>	<i>401mm</i>	<i>401mm</i>
<i>Case DIMENSION (WIDTH)</i>	<i>250mm</i>	<i>250mm</i>	<i>250mm</i>
<i>Case DIMENSION (HEIGHT)</i>	<i>235mm</i>	<i>235mm</i>	<i>235mm</i>
<i>All weights are approximate. make sure to get detailed information from the seller before shipment.</i>			

#### **1.2) Number in dispenser: 100 pcs**

#### **2) product Information:**

- vinyl examination glove
- Powder free
- Excellent finishing
- medical grade
- AQL 1.5



ISO 9001  
ISO 13485  
ISO 10002 & 10004  
FDA QSR 21 CFR 820 (cGMP)  
GDP Good Distribution Practices  
CE (Directive 93/42/EEC, Class I)


- Type of product: Examination Gloves • Material of product: Vinyl (Poly Vinyl Chloride)
- Specification: Non-Sterile, Latex Free & Powder Free, Easy wear, Single use only, Smooth surface
- Number in a box: 100 Pieces by weight • Storage: Store at room temperature, avoid excessive heat (5 - 40° C), Pack should be shield from sunlight, fluorescent lighting, X-ray, moisture and ozone
- Germany Identifier: SONNHEIM Medikal (Germany Reg. No. 39634985.4)
- Manufacturer name & address: Iran Harir Glove Co. / Qazvin Factory: Hafez St., Dehkhoda boulevard Caspian industrial city, Qazvin - I.R.I
- **EC REP TAI Trading Agency GmbH**  
Schüttelstraße 55/11, 1020 Wien-AUSTRIA  
Tel.: +43 (0) 732 77 33 20 / Fax: +43 (0) 732 77 33 20-15  
E-Mail: info@tai-austria.com
- Customer Service (Medical Dep.)  
Tel: + 98 (21) 88 91 64 36 Fax: + 98 (21) 88 91 64 37  
Zip Code: 1415846351 / www.glove.ir / medical@glove.ir  
Made in I.R.I

2.2) CAUTION: -

2.3) symbol used:



EN ISO 374-1: 2016  
type B



LOT



EXP



3) information for transportation:

product name	capacity on pallet (110*130 cm)	capacity in transit trailer, without pallet(approx.)	capacity in transit trailer, with pallet(approx.)
vinyl Examination glove	117	3200	2340

product name	capacity on pallet (110*130 cm)	capacity in 40-foot container, without pallet(approx.)	capacity in 40-foot container, with pallet(approx.)
vinyl Examination glove	104	2430	1872

OP-PERFECT®  
Premium Vinyl Examination Gloves  
Powder Free  
Size M

OP-PERFECT®  
Premium Vinyl Examination Gloves  
Powder Free  
Size M

OP-PERFECT®  
Premium Vinyl Examination Gloves  
Powder Free  
Size M

OP-PERFECT®  
Premium Vinyl Examination Gloves  
Powder Free  
Size M

OP-PERFECT®  
Premium Vinyl Examination Gloves  
Powder Free  
Size M



OP-PERFECT®  
Premium Vinyl Examination Gloves  
Powder Free

Iran Harir Glove Co. Since 1977  
 • Type of product: Examination Gloves • Material of product: Vinyl (Poly Vinyl Chloride)  
 • Specification: Non-sterile, Latex Free & Powder Free, Easy wear, Single use only, Smooth surface  
 • Number in a box: 100 Pieces by weight • Storage: Store at room temperature, avoid excessive heat (5-40° C). Pack should be shield from sunlight, fluorescent lighting, X-ray, moisture and ozone  
 • Germany Identifier: SONNHEIM Medikal (Germany Reg. No. 30634985-4)  
 • Manufacturer name & address: Iran Harir Glove Co. / Gazvin Factory: Hafez St., Dehkhoda boulevard  
 • TAI Trading Agency GmbH  
 Schuttensteinstr. 55/11, 1020 Wien-AUSTRIA  
 Tel.: +43 (0) 732 77 33 20 | Fax: +43 (0) 732 77 33 20-15  
 E-Mail: info@tai-austria.com  
 • Customer Service (Medical Dep.)  
 Tel.: +98 (21) 88 91 64 36 Fax: +98 (21) 88 91 64 37  
 Zip Code: 1415848351 / www.glove.ir / medical@glove.ir  
 Made in Iran



CE  
 ISO 9001  
 ISO 13485  
 ISO 10993 & 10994  
 FDA GSR 21 CFR 820 (cGMP)  
 GMP Good Distribution Practices  
 CE (Directive 93/42/EEC, Class II)



شماره پانک استعمال اصالت پروانه استاندارد ۱۵۷۷-۱۰۰۰  
 کارخانه قزوین: خیابان صنعتی، بلوار دهخدا  
 کارخانه اصفهان: بلوار امام خمینی، خیابان هفتم تیر  
 آدرس: تهران، خیابان هفتم تیر، پلاک ۱۰۰۰  
 شماره تماس: ۰۲۱-۸۸۹۱۶۴۳۶  
 وبسایت: www.glove.ir

000614