

Packing plan for OP-PERFECT EXCELLENT

All issues- warning notices, signs are used, identification plate, storage, transportation, and supply conditions in the plan package are available.

Product Name: latex examination gloves

Brand: OP-PERFECT

1) Packing items:

The initial packaging of the examination gloves is in a Dispenser.

1.1) Packaging information:

<i>NAME</i>	<i>LATEX 100 PCS Size L</i>	<i>LATEX 100 PCS Size M</i>	<i>LATEX 100 PCS Size S</i>
<i>Dispenser WEIGHT (GROSS) (100 pcs)</i>	<i>647 gr</i>	<i>570 gr</i>	<i>554 gr</i>
<i>Case WEIGHT (GROSS) (1000 pcs)</i>	<i>6910 gr</i>	<i>6200 gr</i>	<i>6130 gr</i>
<i>Case DIMENSION (LENGTH)</i>	<i>405mm</i>	<i>405mm</i>	<i>405mm</i>
<i>Case DIMENSION (WIDTH)</i>	<i>255 mm</i>	<i>255 mm</i>	<i>255 mm</i>
<i>Case DIMENSION (HEIGHT)</i>	<i>235mm</i>	<i>235mm</i>	<i>235mm</i>
<i>All weights are approximate. make sure to get detailed information from the seller before shipment.</i>			

1.2) Number in Dispenser: 100 pcs

2) product Information:

- latex examination glove*
- Powder free*
- Excellent finishing*
- medical grade*
- AQL 1.5*



Iran Harir Glove Co.
Since 1977



ISO 9001
ISO 13485
ISO 10002 & 10004
FDA QSR 21 CFR 820 (cGMP)
GDP Good Distribution Practices
CE (Directive 93/42/EEC, Class I)

- Type of product: **Examination Gloves** ● Material of product: **Latex (NR)**
- Specification: **Non-Sterile, Latex & Powder Free, Easy wear, Single chlorinated, Polymer plus, Ambidextrous Micro textured, Single use only** ● Number in a box: **100 Pieces by weight**
- ⚠ Cautions: **This product contains natural rubber latex which may cause allergic reactions in some individuals**
- Storage: **Store at room temperature, avoid excessive heat (5- 40°C), Pack should be shield from sunlight, fluorescent lighting, X-ray, moisture and ozone**
- Germany Identifier: **SONNHEIM Medikal (Germany Reg. No. 39634985.4)**
- 🏭 Manufacturer name & address: **Iran Harir Glove Co. Urmia Factory: Phase #2 of Industrial city, Urmia - I.R.I**
- EC REP TAI Trading Agency GmbH / Schüttelstraße 55/11, 1020 Wien-AUSTRIA / Tel.: +43 (0) 732 77 33 20 Fax: +43 (0) 732 77 33 20-15 / E-Mail: info@tai-austria.com**
- Customer Service (Medical Dep.) Tel: **+ 98 (21) 88 91 64 36 Fax: + 98 (21) 88 91 64 37**
- www.glove.ir / medical@glove.ir / Zip Code: 1415846351 Made in I.R.I**

2.2) CAUTION:

⚠ Cautions: **This product contains natural rubber latex which may cause allergic reactions in some individuals**

2.3) symbol used:



LOT



EXP



3) information for transport:

Product name	Capacity on pallet (110 * 100 cm)	Capacity in transit trailer, without pallet (approx.)	Capacity in transit trailer, with pallet (approx.)	Capacity in 40-foot container, without pallet (approx.)	Capacity in 40-foot container, with pallet (approx.)
powder free Examination glove	90	3260	2340	2860	1980

OP-PERFECT Powder Free Examination Glove
EXCELLENT FINISHING
Size M

Latex
Powder Free

100 Pieces by Weight

CE ISO 9001 ISO 13485 ISO 14685 & 14686 FDA QSR 21 CFR 820 (cGMP) GMP Good Distribution Practices CE (Directive 93/42/EEC, Class II)

OP-PERFECT Powder Free Examination Glove
EXCELLENT FINISHING
Size M

OP-PERFECT Powder Free Examination Glove
EXCELLENT FINISHING
Size M

OP-PERFECT Powder Free Examination Glove
EXCELLENT FINISHING
Size M

Iran Harir Glove Co. Since 1977

6 260045 229949

Size M
Medical Grade
AQL 1.5

XS S M L XL XXL

XS S M L XL XXL

• Type of product: Examination Glove • Material of product: Latex (NR)
 • Specification: Non-Sterile, Latex & Powder Free, Easy wear, Single chlorinated, Polymer plus, Antistatic
 Micro textured, Single use only • Number in a box: 100/Pieces by weight
 A. Caution: This product contains natural rubber latex which may cause allergic reactions in some individuals
 • Storage: Store at room temperature, avoid excessive heat (5-40°C). Pack should be shield from sunlight, fluorescent lighting, X-ray, moisture and ozone
 • Manufacturer name & address: Iran Harir Glove Co
 Urmia Factory: Phase #2 of Industrial city, Urmia - I.R.
 [EU] 2828 1st Trading Agency GmbH / Schindelsiedle 65/11, 1020 Wien-AUSTRIA / Tel.: +43 (0) 732 77 33 20
 Fax: +43 (0) 732 77 33 20-12 / E-Mail: info@rai-austria.com
 • Customer Service (Medical Dept) Tel.: +98 (0) 88 91 64 38 Fax: +98 (0) 88 91 64 37
 www.glove.ir / medical@glove.ir / Zip Code: 1419848351 Made in I.R.