

Specification for OP-PERFECT vinyl examination gloves

color	<i>NATURAL COLOR (Multicolor on request)</i>								
shape	Smooth								
cuff	Beaded								
Polyurethane Coated Powder Free Vinyl Examination Gloves									
PHYSICAL DIMENSIONS									
Size	Length (mm)		Palm Width (mm)		Thickness(mm)				
	Standard	Test	Standard	Test	Standard	Test			
						Cuff	Palm	Fin	
Small	≥230	248	87±3	87	Textured area 0.11 mm	0.08	0.08	0.11	
Medium		246	99±5	99		0.08	0.08	0.12	
Large		241	110±5	106	Smooth area 0.08 mm	0.08	0.08	0.12	
Extra Large		243	≥110	116		0.08	0.08	0.10	
Tensile Properties	Standard		Test						
			1	2	3	4	5	6	7
Before aging	Tensile(N)	7.0	7.31	7.26	7.42	7.19	7.37	7.52	7.23
	Elongation (%)	350	385	379	396	383	391	395	379
After aging	Tensile (N)	7.0	7.06	7.13	7.21	7.03	7.17	7.28	7.0
	Elongation (%)	350	357	366	382	378	371	383	362
Protein(μgr/gr)		Not applicable							
Powder(mg)		Not applicable							