

Op-perfect HARIR SURGICAL GLOVES

Op-perfect HARIR

Latex, powder-free surgical glove offering unsurpassed combination of comfort, protection and sensitivity

Driven by innovation, OP-PERFECT surgical gloves are engineered to respond to the evolving users needs of the surgical community. Featuring unique formulations and cutting-edge technologies, ergonomical comfort and higher-level protection, healthcare workers around the world entrust their hands to us.

KEY FEATURES AND BENEFITS:

- Optimized glove thickness for the finest balance of protection and tactile sensitivity*
- Additional cleaning processes to minimize the risk of developing allergies and sensitivities*
- Optimal tactile sensitivity and protection*

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

Op-perfect HARIR SURGICAL GLOVES

PRODUCT INFORMATION

<i>Material</i>	<i>Natural Rubber Latex</i>	<i>Not Made from Natural Rubber Latex</i>	<i>No</i>
<i>Color</i>	<i>Natural color</i>	<i>Cuff Style</i>	<i>Beaded</i>
<i>Allergy Prevention</i>	<i>None</i>	<i>Sterile</i>	<i>yes</i>
<i>Powder Content</i>	<i>POWDER-FREE</i>	<i>outer Surface</i>	<i>micro-Textured</i>
<i>Freedom from Holes (Inspection level I)</i>	<i>0.65 AQL</i>	<i>Palm Thickness (mm)</i>	<i>0.16-017</i>
<i>Finger Thickness (mm)</i>	<i>0.22-0.23</i>	<i>Available Sizes</i>	<i>6.5-7-7.5-8-8.5</i>
<i>Inner Surface</i>	<i>Polymer Coated</i>	<i>Glove Length (mm)</i>	<i>285-305</i>
<i>Product Standards</i>	<i>EN 455 1-4, ISO 10282</i>	<i>Quality / Environmental Standards</i>	<i>EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001</i>

ORDERING INFORMATION

<i>Size</i>	<i>6.5</i>	<i>7</i>	<i>7.5</i>	<i>8</i>	<i>8.5</i>
<i>Product Code</i>	<i>741011132002</i>	<i>741011132003</i>	<i>741011132004</i>	<i>741011132005</i>	<i>741011132006</i>

PACKAGING AND STORAGE

<i>Packaging</i>	<i>400 pairs per case; 8 dispensers per case; 50 pairs/envelopes per dispenser</i>
<i>Shelf Life</i>	<i>5 Years</i>
<i>Storage Instructions</i>	<i>Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition</i>