

Packing plan for OP-PERFECT PREMIUM

All issues- warning notices, signs are used, identification plate, storage, transportation, and supply conditions in the plan package are available.

Product Name: latex surgical gloves

Brand: OP-PERFECT PREMIUM

1) Packing items:

The initial packaging of the surgical gloves is in an envelope.

1.1) Packaging information:

<i>NAME</i>	<i>LATEX 50 pairs Size 6.5</i>	<i>LATEX 50 pairs Size 7</i>	<i>LATEX 50 pairs Size 7.5</i>	<i>LATEX 50 pairs Size 8</i>	<i>LATEX 50 pairs Size 8.5</i>
<i>Dispenser WEIGHT (GROSS) 50 pairs</i>	<i>1740 gr</i>	<i>1780gr</i>	<i>1850 gr</i>	<i>1960 gr</i>	
<i>Case WEIGHT (GROSS) 400 pairs</i>	<i>14670 gr</i>	<i>15000gr</i>	<i>15600 gr</i>	<i>16480 gr</i>	
<i>Case DIMENSION (LENGTH)</i>	<i>620mm</i>	<i>620mm</i>	<i>620mm</i>	<i>620mm</i>	
<i>Case DIMENSION (WIDTH)</i>	<i>296 mm</i>	<i>296 mm</i>	<i>296 mm</i>	<i>296 mm</i>	
<i>Case DIMENSION (HEIGHT)</i>	<i>400mm</i>	<i>400mm</i>	<i>400mm</i>	<i>400mm</i>	
<i>All weights are approximate. make sure to get detailed information from the seller before shipment.</i>					

1.2) Number in envelope: 1 pair

2) product Information:

- latex surgical glove*
- low powder*
- Excellent finishing*
- medical grade*

- **Type of product:** Surgical Gloves
 - **Material of product:** Latex (NR)
 - **Specification:** Sterile, Latex & Low Powder, Easy wear, Single chlorinated, Single use only • Palm textured
 - The **protein** level of this product is lower than ASTM D 5712-15 Standard
 - Every glove 100% is **air inflation** tested and **visually** inspected for quality and safety
 - Improved **efficiency** as less glove are wasted
 - Recommended for **all surgical procedures**
 - **Extra Length** for more protection
 - Disposable, Soft, Sensible
 - **Anatomic** shape
 - **Non-beaded** cuff
 - Improved **grip** to prevent slipping **surgical instrument**
 - **Sterilized by:** Ethylene oxide (E.O) gas
 - Number in a pack: **1 Pair**
- ⚠ **Cautions:** This product contains natural rubber latex which may cause allergic reactions in some individuals
- Surface powder should be removed after donning
 - ☯ Sterile unless package is opened or damaged
- **Storage:** Store at room temperature, avoid excessive heat (5 - 40°C), Pack should be shield from sunlight, fluorescent lighting, x-ray, moisture and ozone
- **Identifiers:** ISO 13485 / ISO 9001 / ISO 10002 & 10004 / FDA QSR 21 CFR 820 (cGMP)
GDP Good Distribution Practices / ISIRI (Cer. No. 6117041877)
Health sign of Ministry of Health (Cer. No. 68937701) / SONNHEIM Medikal (Germany Reg. No. 39634985.4)
- 🏭 **Manufacturer name & address:** Manufacturer name & address: Iran Harir Glove Co. / Urmia Factory: Phase #2 of Industrial city, Urmia - I.R.I
- ECIREP TAI Trading Agency GmbH / Sch üttelstra e 55/11, 1020 Wien-AUSTRIA / Tel.: +43 (0) 732 77 33 20
Fax: +43 (0) 732 77 33 20-15 / E.Mail: info@tai-austria.com
- **Customer Service (Medical Dep.)** Tel: + 98 (21) 88 91 64 36 Fax: + 98 (21) 88 91 64 37
www.glove.ir / medical@glove.ir / Zip Code: 1415846351 / Made in I.R.I

2.2) CAUTION:

- ⚠ **Cautions:** This product contains natural rubber latex which may cause allergic reactions in some individuals
- Surface powder should be removed after donning
 - ☯ Sterile unless package is opened or damaged

2.3) symbol used:



1Pair
STERILE EO

Size

LOT



3) *information for transportation:*

<i>Product name</i>	<i>Capacity on pallet (110 * 130 cm)</i>	<i>Capacity in transit trailer, without pallet (approx.)</i>	<i>Capacity in transit trailer, with pallet (approx.)</i>	<i>Capacity in 40-foot container, without pallet (approx.)</i>	<i>Capacity in: 40-foot container, with pallet (approx.)</i>
<i>Latex surgical gloves</i>	<i>35</i>	<i>1050</i>	<i>770</i>	<i>756</i>	<i>630</i>



Latex, Low Powder Surgical Gloves **OP-PERFECT**[®] PREMIUM FINISHING



Iran Harir Glove Co.
Since 1977



- **Type of product:** Surgical Gloves
- **Material of product:** Latex (NR)
- **Specification:** Sterile, Latex & Low Powder, Easy wear, Single chlorinated, Single use only • Palm textured
- The **protein** level of this product is lower than ASTM D 5712-15 Standard
- Every glove 100% is **air inflation** tested and **visually** inspected for quality and safety
- Improved **efficiency** as less glove are wasted
- Recommended for **all surgical procedures**
- **Extra Length** for more protection
- Disposable, Soft, Sensible
- **Anatomic** shape
- **Non-beaded** cuff
- improved **grip** to prevent slipping surgical instrument
- **Sterilized by:** Ethylene oxide (E.O) gas
- Number in a pack: **1 Pair**

⚠ **Cautions:** This product contains natural rubber latex which may cause allergic reactions in some individuals

- Surface powder should be removed after donning
- Sterile unless package is opened or damaged

• **Storage:** Store at room temperature, avoid excessive heat (5 - 40°C), Pack should be shield from sunlight, fluorescent lighting, x-ray, moisture and ozone

• **Identifiers:** ISO 13485 / ISO 9001 / ISO 10002 & 10004 / FDA QSR 21 CFR 820 (cGMP)
GDP Good Distribution Practices / ISIRI (Cer. No. 6117041877)

Health sign of Ministry of Health (Cer. No. 68937701) / SONNHEIM Medikal (Germany Reg. No. 39634985.4)
• **Manufacturer name & address:** Manufacturer name & address: Iran Harir Glove Co. / Urmia Factory: Phase #2 of Industrial city, Urmia - I.R.I

• **REGD:** TAI Trading Agency GmbH / Sch ü ttelstra e 55/11, 1020 Wien-AUSTRIA / Tel.: +43 (0) 732 77 33 20
Fax: +43 (0) 732 77 33 20-15 / E.Mail: info@tai-austria.com

• **Customer Service (Medical Dep.)** Tel: + 98 (21) 88 91 64 36 Fax: + 98 (21) 88 91 64 37
www.glove.ir / medical@glove.ir / Zip Code: 1415846351 / **Made in I.R.I**

تولید کننده: شرکت صنایع دستکش ایران
تولید کننده: شهر صنعتی ارومیه، شهر صنعتی، فاز ۲
تولید کننده: شهر صنعتی ارومیه، خیابان
تولید کننده: شهر صنعتی ارومیه، خیابان
تولید کننده: شهر صنعتی ارومیه، خیابان
تولید کننده: شهر صنعتی ارومیه، خیابان



شماره پیامک استاندارد پروانه استاندارد ۱۰۰۰۱۵۱۷

معاف از مالیات بر ارزش افزوده

Size

LOT

